

TEST PROCEDURE

Check instructions Use for complete procedue

1



2




3



4



5



CONTACT US

BC MEDICARE SDN BHD

Level 3, Tower 11, Avenue 5, No.8 Jalan Kerinchi
Bangsar South. 59200 Kuala Lumpur, Malaysia

Website:

<https://www.bcmalliance.com.my>

Email:

sghoo@bcmalliance.com.my

Tel: +603-2703 3753

Whatsapp: +6010- 393 4135



SARS-COV-2 Antigen Rapid Test Kit

less-invasive COVID-19 testing
with nasal (NS) swab

BC MEDICARE SDN BHD
Website: <https://www.bcmalliance.com.y>

PRODUCT INFO

* Performance:

- SENSITIVITY: 98.13%
- SPECIFICITY: 99.22%

SPECIFICATION:

TEST TIME: 15-20 MINUTES
STORAGE CONDITON: 2°C- 30°C
CE MARK
SAMPLE TYPE: NASAL (NS) SWAB
SHELF LIFE: 24 MONTHS

CONTENTS:

20x Test Devices
20x Buffer Bottle (350µl/bottle)
20x Extraction Tubes & Nozzle Caps
1x Tube Rack
1x Instruction for Use
Positive/Negative Control Swabs
(upon request)

HIGHLIGHT

High performance proved by an independent clinical study in Europe
Rapid detection of SARS-COV-2 nucleocapsid antigen in less than 15 minutes

Less-invasive sample collection with nasal swabs.

No instrumentation required
Cost-effective solution for large-scale testing.



INTENDED USE

For in vitro qualitative detection of SARS-CoV-2 nucleocapsid antigen in nasal (NS) swabspecimens directly from individuals who meetCOVID-19 clinical and / or epidemiologicalcriteria. SARS-CoV-2 Antigen Rapid TestKit is for professional use only and isintended to be used as an aid in the diagnosisof SARS-CoV-2 infection. The product may beused in any laboratory and non-laboratoryenvironment that meets the requirementsspecified in the Instructions for Use and local regulation.



* Sensitivity and specificity vs PCR using Nasopharyngeal swabs.